## **Kable EKG do badań holterowskich**

**ECG holter cables**

**Langzeit - EKG Kabel**

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EN

## Instruction of use and storage ECG holter cables

The instruction is concerned with ECG holter cables for type KB, KH, KM, KS

### 1. General description

Cables are intended to transmit electrical signals received from human organism to devices processing these signals - ECG apparatus.

Cables are non-sterile class I medical products.

In order to ensure proper cooperation of ECG cable with a monitor, choose an appropriate cable compatible with the device.

### 2. ECG cable characteristics

The amount of conductors in the cable, cable length, method of connection of individual conductors to connectors, type of connectors, placement of electronic elements such as resistors and anti-surge elements in the cable, are adjusted to a specific type/model of apparatus, with which the cable is to cooperate.

Marking of the cable, method of its connection, purpose of individual tips and outlets, and other data, essential for correct conduction of examinations, is specified in the instruction of use of the device, with which the cable is intended to cooperate.

Cables are intended for use with undamaged patient's skin.

Cables should be used exclusively by qualified personnel.

### 3. Principles of use

When using ECG cables, the following principles should be observed:

- cables should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the cable,
- place the cable so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin, clasp, latch ) when connecting and disconnecting the cable,
- protect the cable connectors against suffusion.

### 4. Installation

Connect the cable pin to the appropriate ECG apparatus, turn on the device and verify the correctness of functioning according to monitor's instruction of use.

Pay attention to proper choice of the cable for an appropriate ECG apparatus.

### 5. Recommended cleaning, disinfection, and sterilisation

These activities should be realized every time before use of the device and cable before the examination of the next patient. Before starting the cleaning, disinfection, or sterilisation, the cable should be disconnected from the device with which it cooperates.

#### Cleaning

Wash the cable with warm soap water, delicately clean using soft cloth wetted with soap water

(mild detergent) or 70% isopropyl alcohol. The cable should dry before use.

Do not use corrosive substances or other chemical preparations.

When cleaning, do not soak or immerse the cables, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to cable damage. Continuous bending of the conductors during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.

Avoid contact with strong solvents: aromatic, chloric, ketone, etheral, and ester. Do not allow for contact of cleaners with tips from the electrode side, and pins from ECG apparatus side.

#### Disinfection

Cables should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry. Do not allow for contact of alcohol with tips from the electrode side, and pins from ECG apparatus side.

#### Sterilization

The recommended method consists of sterilization using ethylene oxide. After sterilization and before use, cables should be accurately degassed in device intended for that purpose.

Sterilization with other methods is not recommended for it may permanently damage the cable.

### **6. Packaging and storing**

ECG cables are packed individually. In order to maximize the product's working life, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:

- temperature: -25°C to +55°C,
- relative humidity: 10% to 85%, do not allow for steam condensation
- protect against sunlight, moisture and dirt.

Cables in use should be stored in a dry place.

When not used, ECG cables should be loosely coiled, without sharp bends.

### **7. Environment conditions during work**

- temperature +5°C to +45°C
- relative humidity 30% to 60%
- atmospheric pressure 700 ÷ 1060 hPa

### **8. Warnings**

Not observing the principles provided above may break electric circuit or lead to short-circuit. These signs are easy to identify by the ECG apparatus service staff. Damaged cable should be put away and replaced with an operational one.

ECG cables are intended for use with specified types of ECG apparatuses. The user is responsible for verification, before its use, of cable compatibility with the device. Incompatibility of these elements may lead to reduced accuracy and efficiency of measurements (or lack of them), and in extreme case - to damaging of the product. In order to ensure their compatibility, check the equipment cooperating with a given device in its instruction of use.

It is not recommended to use ECG cable in aroused patients, for movement negatively influences realized measurements. Relocation and positioning of the patient may cause disconnection of the cable from ECG electrodes and lack of readings.

Strong sources of electromagnetic field, such as those emitted by electrosurgical equipment, may negatively influence the correct functioning of the ECG cable.

Allergic reactions at contact of ECG cable with patient's undamaged skin are not observed. However, such reactions may not be excluded for persons particularly prone to allergies. Improper cleaning or disinfection (or lack of them), and disregard for hygienic procedures may result in transfer of micro-organisms and cause infections.

Do not use an ECG cable suspected of having mechanical or electrical defects, or packaging of which has been damaged. Suffused ECG cable should not be used. Do not change or modify the cable in any way. Such actions may influence correctness of functioning and accuracy of

received readings.

Utilisation of used cable should be done in accordance with binding provisions.

## 9. Guarantee

The guarantee period covers 1 year from the date of sale.

Shelf life is 2 years from the date of sale.

The guarantee does not include defects resulting from improper use and storage.

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