SORIMEX

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Leadwires

Verbindungskabel

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EN

Instruction of use and storage Leadwires

The instruction is concerned with leadwires type PD and PG

1. General description

Leadwires are intended to transmit electrical signals received from human organism to devices processing these signals - ECG apparatus, unassisted or with consolidated cable.

Leadwires are non-sterile class I medical products.

2. Leadwires characteristics

Type of applied wire, its length, type of endings on patient's side, usage of resistance or inductive elements are adjusted to specific type/model of apparatus or consolidated cable with which the leadwire is to cooperate.

Marking of the leadwire, method of its connection, purpose of individual tips and outlets, and other data essential for correct conduction of examinations, is specified in the instruction of use of the device to cooperate with which the leadwire is intended.

Leadwires are intended for use with undamaged patient's skin.

Leadwires should be used exclusively by qualified medical personnel.

3. Principles of use

When using leadwires, the following principles should be observed:

- leadwires should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the leadwire,
- place the leadwire so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin, clasp, latch) when connecting and disconnecting the leadwire,
- protect the leadwire connectors against suffusion.

4. Installation

Connect the cable pin to the appropriate consolidated cable/ECG apparatus, turn on the device and verify the correctness of functioning according to monitor's instruction of use.

Pay attention to proper choice of the leadwire for appropriate consolidated cable/ECG apparatus.

5. Recommended cleaning, disinfection, and sterilisation

These activities should be realized every time before use of the device and leadwire before the examination of the next patient. Before starting the cleaning, disinfection, or sterilisation, the leadwire should be disconnected from the device with which it cooperates.

Cleaning

Wash the leadwire by gently wiping it with soft cloth moistened with warm water with soap or mild detergent; the leadwire should be left to dry before use.

Do not use corrosive substances or other chemical preparations.

When cleaning, do not soak or immerse the leadwires, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to leadwire

damage. Continuous bending of the conductors during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.
Avoid contact with strong solvents: aromatic, chloric, ketone, etheral, and ester. Do not allow for contact of cleaners with tips from the electrode side, and pins from ECG apparatus side.

Disinfection

Leadwires should be disinfected by wiping with 70% isopropyl alcohol and leaving them to dry. Do not allow for contact of alcohol with tips from the electrode side, and pins from ECG apparatus side.

Sterilization

The recommended method consists of sterilization using ethylene oxide. After sterilization and before use, leadwires should be accurately degassed in device intended for that purpose. Sterilization with other methods is not recommended for it may permanently damage the leadwire.

6. Packaging and storing

Leadwires are packed individually. In order to maximize the product's working life, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:

- temperature: -25°C to + 55°C,
- relative humidity: 10% to 85%, do not allow for steam condensation
- protect against sunlight, moisture and dirt.

Leadwires in use should be stored in a dry place.

When not used, leadwires should be loosely coiled, without sharp bends.

7. Environment conditions during work

- temperature +5°C to +45°C
- relative humidity 30% to 60%
- atmospheric pressure 700 + 1060 hPa

8. Warnings ⚠

Not observing the principles provided above may break electric circuit or/and leadwire shield, or lead to short-circuit. These signs are easy to identify by the ECG apparatus service staff. Damaged leadwire should be put away and replaced with an operational one.

Leadwires are intended for use with specified types of ECG apparatuses. The user is responsible for verification, before its use, of leadwire compatibility with the device. Incompatibility of these elements may lead to reduced accuracy and efficiency of measurements (or lack of them), and in extreme case - to damaging of the product. In order to ensure their compatibility, check the equipment cooperating with a given device in its instruction of use.

It is not recommended to use leadwire in aroused patients, for movement negatively influences realized measurements. Relocation and positioning of the patient may cause disconnection of the leadwire from ECG electrodes and lack of readings.

Strong sources of electromagnetic field, such as those emitted by electrosurgical equipment, may negatively influence the correct functioning of the leadwire.

Allergic reactions at contact of leadwire with patient's undamaged skin are not observed. However, such reactions may not be excluded for persons particularly prone to allergies. Improper cleaning or disinfection (or lack of them), and disregard for hygienic procedures may result in transfer of micro-organisms and cause infections.

Do not use an leadwire suspected of having mechanical or electrical defects, or packaging of which has been damaged. Suffused leadwire should not be used.

Do not change or modify the leadwire in any way. Such actions may influence correctness of functioning and accuracy of received readings.

Utilisation of used leadwire should be done in accordance with binding provisions.

9. Garantie

The guarantee period covers 1 year from the date of sale.

Shelf life is 2 years from the date of sale.

The guarantee does not include defects resulting from improper use and storage.

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