# SORIMEX

## **PL Elektrody EKG w.u.**

## **EN Reusable ECG electrodes**

SORIMEX sp. z o.o. (dawniej SORIMEX sp. z o.o. sp. k.)  
Równinna 25, 87-100 Toruń, Poland  
office@sorimex.eu ♦ www.sorimex.eu

Rev. 2025-01-07

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PL

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EN

This instruction for use concerns:

- **Reusable ECG electrodes** of the following types:
  - EKK
  - EP

## 1. Intended use and operation of the product

The electrodes are intended to transmit the electrical signals from the body surface through ECG leads and cables to the ECG unit or a cardiac monitor, which processes the obtained signal into an electrocardiogram or a vectorcardiogram.

## 2. Target patient groups

Electrodes are intended to be used by qualified medical personnel in children and adults.

## 3. Contraindications

No known contraindications to the use of the product.

## 4. Possible adverse effects

Electrodes are devices having a direct contact with the patient's skin. There is a risk of irritation and allergy reactions in patients with particularly sensitive skin. If any skin reactions develop, medical advice should be sought.

## 5. Warnings

- 5.1. Disconnect a damaged electrode and replace it with one in a good operating condition. EKK electrodes must be used with units that are in a perfect operating condition.
- 5.2. Apply the electrodes solely to the uninjured skin. Do not apply the electrodes on the skin that is injured, infected or covered with lesions.
- 5.3. Do not use electrodes during magnetic resonance imaging (MRI) scans.
- 5.4. Do not open the electrode outside the extent resulting from its design, so the spring component of the electrode does not break.
- 5.5. To connect the electrode with the ECG unit use cables of a design suitable for operation with a given ECG unit.

## 6. Safety precautions

- 6.1. Do not use electrodes that have mechanical defects.
- 6.2. Do not allow the electrode to contact the injured skin.
- 6.3. Excessive movements of the patient adversely affect the quality of the ECG signal.

## 7. How to use ECG electrodes

- 7.1. Determine places for attachment of electrodes on the patient's body

- 7.2. Clean the skin in selected area of creams, ointments, etc. If necessary, degrease the patient's skin with water and soap, and dry thoroughly.
- 7.3. Apply gel for ECG examinations on the patient's skin in places in contact with an electroconductive part.
- 7.4. Connect to leads of the ECG cable.
- 7.5. Check the electrode - cable - unit connection for its durability.
- 7.6. Ensure that the ECG cable is placed in a way reducing a possibility of patient's entangling or of damaging the product.

## **8. Cleaning and disinfection**

These actions should be performed each time before the electrodes are used.

### **8.1. Cleaning**

Electrodes should be cleaned by delicately wiping them with a soft cloth wetted with warm water and soap or mild detergent. The electrode should be dry before use.

During cleaning, do not wet or immerse the electrode, do not wash under running water, and do not use excessive force that may damage the electrode.

Do not use any corrosive or abrasive agents or other aggressive chemical formulations.

Avoid any contact with strong solvents: aromatic, chlorine, ketone, ethers or esters.

### **8.2. Disinfection**

Disinfect the electrodes by wiping with formulations intended for disinfection of medical devices made of plastic, if they do not negatively affect the electrode. Be sure to read the instructions for use of those formulations carefully before their use.

Formulations recommended for disinfection of the electrodes include 70% of isopropyl alcohol or 70% ethanol.

## **9. Signs of product degradation**

Before each use, perform visual inspection of the product. Electrodes cannot be used any longer when any of the below is noticed:

- interrupted continuity of the plastic surface;
- damaged or corroded sensors.

Any of the above defects disqualifies the product from further use. A product with such defects should either be repaired or disposed of.

## **10. Storage conditions**

### **10.1. Storage conditions:**

- temperature: 5°C to +30°C
- protect against sunlight, moisture and contamination

**10.2.** Electrodes should be stored in the specified storage conditions, to prolong the product life as far as possible.

## 11. Compatibility

Electrodes are connected to the ECG unit with an ECG cable ending with a banana plug of  $\varnothing$  4 mm.

## 12. Warranty

The product life is 5 years of the production date specified in the label.

## 13. Description of symbols used

Product catalogue number

Pieces per packaging

Production batch number

Medical device

Production date

Manufacturer

Consult the instruction for use

Caution

Protect against sunlight

CE mark  
(conformance to the Regulation EU 2017/745)

Keep away from rain

Acceptable temperature limits

Product index

Product UDI code

## 14. Notification of incidents related to the product use

Each serious incident related to the product should be notified to the manufacturer and to an authority of the Member State competent for the user or patient's place of residence.

## 15. Product disposal

The electrodes should be disposed of in accordance with internal hygiene procedures in force at a given health care facility and with local regulations.