# SORIMEX

## **PL Kable czujników temperatury**

## **EN Temperature sensor cables**

SORIMEX sp. z o.o. sp. k.  
Równinna 25, 87-100 Toruń, Poland  
office@sorimex.eu ♦ www.sorimex.eu

Rev. 2022-04-06

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PL

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- **kabli czujników temperatury**

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This instruction for use concerns:

- **Temperature sensor cables**

## 1. Intended use and operation of the product

The temperature sensor cable is intended for electrical connection of the temperature sensor with the device.

The cable consists of a multi-wire cable of an appropriate length ending with connectors suitable for the temperature sensor and the device with which it is to work.

## 2. Target patient groups

Temperature sensor cables are intended to be used by qualified medical personnel in children and adults.

## 3. Contraindications

No known contraindications to the use of the product.

## 4. Possible adverse effects

Temperature sensor cables can contact the uninjured patient's skin. No allergic reactions are observed in contact with the patient's skin; however, such reactions cannot be completely excluded in patients particularly susceptible to allergies.

## 5. Warnings

- 5.1. A damaged temperature sensor cable should be disconnected and replaced with one in a good operating condition. Temperature sensor cables must be used with units that are in a perfect operating condition.
- 5.2. Temperature sensor cables are intended to be used with specific types of temperature sensors and devices. Before each use of the cable, check whether it is suitable for a given type of the device and the sensor, with which it is to operate.
- 5.3. When installing the temperature sensor cable, follow the operating manual of the device.

## 6. Safety precautions

- 6.1. Do not use temperature sensor cables with any mechanical or electrical defects.
- 6.2. Do not use temperature sensor cables that were in contact with any liquid.
- 6.3. Temperature sensor cables should not be bent excessively (the bend radius should be at least 60 mm).
- 6.4. Do not pull or yank a cable.
- 6.5. A cable should be routed in such way that there is no risk of its mechanical damage caused, e.g., by stepping on, driving over, placing a heavy/sharp object on, or cutting it.

- 6.6. When connecting or disconnecting a cable, hold it by a connection (socket, plug) housing.
- 6.7. Movement, transfer and positioning of the patient may interfere with reception of collected signals and affect measurement results.
- 6.8. Do not use the temperature sensor cable during Magnetic Resonance Imaging (MRI) or Computed Tomography (CT) scans.

## **7. How to use the cable**

- 7.1. Each time before using the product, clean and disinfect it.
- 7.2. Select a temperature sensor cable appropriate for the sensor and the device with which it is to operate.
- 7.3. Connect the temperature sensor cable plug to the device socket, and on the other side, connect the temperature sensor cable socket to a plug of the temperature sensor.
- 7.4. Switch the device on and check it for correct operation in accordance with the device operating manual.

## **8. Cleaning and disinfection**

These actions should be performed each time before the temperature sensor cable is used. Before cleaning or disinfection, a cable should be disconnected from a device and a sensor with which it operates.

### **8.1. Cleaning**

A cable should be cleaned by delicately wiping it with a soft cloth wetted with warm water and soap or mild detergent. Before the use, the cable should be dry.

Do not use any corrosive agents or other aggressive chemical formulations. Avoid any contact with strong solvents: aromatic, chlorine, ketone, ethers or esters.

Do not allow any contact of cleaning agents with a terminal of a cable on the sensor side and with plugs on the device side.

### **8.2. Disinfection**

Disinfect cables by wiping them delicately with formulations dedicated to disinfect plastic medical devices, based on active ingredients such as (70%) isopropyl alcohol, glutaraldehyde, amines, and quaternary ammonium compounds. Disinfect according to instructions of a formulation manufacturer, and leave the cable to dry.

## **9. Signs of product degradation**

Before each use, perform visual inspection of the product. Temperature sensor cables cannot be used any longer when any of the below is noticed:

- breaks in the cable jacket (external insulation) continuity, revealing conducting elements of the cable (wires, shield, connecting layer)
- deep scratches not affecting the jacket continuity, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection
- cracks in plastic material covering product components, making electrical connections visible

- deep scratches in plastic material covering product components that do not make electrical connections visible, but which may cause accumulation of pathogens and prevent effective cleaning and disinfection of the product
- partial or complete tearing of cables from connections/plugs or breaking of the jacket in those areas

Any of the above defects disqualifies the product from further use. A product with such defects should either be repaired or disposed of.

## **10. Storage conditions**

### **10.1. Storage conditions:**

- temperature: -25°C to +55°C
- protect against sunlight, moisture and contamination.

**10.2.** Temperature sensor cables should be stored in the specified storage conditions, to prolong the product life as far as possible.

**10.3.** During storage, the temperature sensor cables should be loosely wound, without kinks.

## **11. Compatibility**

The temperature sensor cables are intended to be used as intermediate components connecting temperature sensors with the devices.

A device operating together with the product is specified in a label on the product primary packaging.

## **12. Warranty**

The product life is 36 years of the production date specified in the label. The guarantee period is 12 months of the sales date.

The general guarantee terms and conditions are available at the manufacturer website, [www.sorimex.pl](http://www.sorimex.pl). The guarantee does not cover defects resulting from a failure to observe provisions of this operating manual.

### 13. Description of symbols used

Catalogue number

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Batch code

Medical device

Date of manufacture

Manufacturer

Consult instructions for use or consult electronic instructions for use

Caution

Keep away from sunlight

CE mark (conformance to the Regulation EU 2017/745)

Keep dry

Temperature limitation

Unique Device Identifier

Marking of electrical and electronic devices according to EN 50419

### 14. Notification of incidents related to the product use

Each serious incident related to the product should be notified to the manufacturer and to an authority of the Member State competent for the user or patient's place of residence.

### 15. Product disposal

The used temperature sensor cable should be disposed of in accordance with internal hygiene procedures in force at a given health care facility and with local regulations.